Example 13.5: B₁₈H₂₂

Last update 12.3.2021

Data	Торіс	Level
Electron diffraction	Complete procedure from data reduction to	Intermediate
	structure refinement including microscope	
	distortions using calibrated values.	

Introduction

In this example, one of the possible ways for obtaining the precise lattice parameters is shown. Because the B₁₈H₂₂ is a molecular crystal it is necessary to not only to count with the distortions of the diffraction data by the microscope optics but also with the beam damage induced changes of the lattice. Molecular crystals are significantly beam senzitive, which results in lattice parameters changes induced and intensified by the accumulated electron dose. The optical distortions induced by the microscope optics were calibrated using Lutetium aluminium garnet. These calibrations can then be used for correction the optical distortions in the data obtained on the molecular crystal materials. The stability of the distortions is relatively good so it can be done only once in a while. The problem of lattice parameters evolution stemming from the accumulated electron dose is a much less defined problem and can be mitigated by using symmetry restrictions following the crystal system. The unit cell refinement using these symmetry restrictions is implemented in PETS 2.0.

The example will also show you difference in kinematical and dynamical refinement. You will see that it is significantly easier to locate hydrogen atoms and the temperature factors are more reflecting the real situation in the material.

Input data

Data

Electron diffraction data were measured on a transmission electron microscope FEI Tecnai G^2 20 (LaB₆) with a precession device (DigiStar). Accelerating voltage was 200 kV. The sample temperature was 100 K.

Input files

folder "dp-100": contains 64 measured diffraction patterns (16-bit TIF files) B18H22.pts: input file for PETS Calibrated_optical_distortions.txt: values of pincushion, spiral and elliptical distortions of the microscope and setting of the distortion refinement keys.

Additional information

Supplementary files for easier reproducibility: folder "PETS_output_CIF": expected output of the data reduction folder "JANA": intermediate states of the refinements. Thus, different steps can be evaluated without having to carry out the computationally expensive (i.e. time consuming) dynamical refinement steps Software used: Jana2020 (jana.fzu.cz) PETS (pets.fzu.cz) VESTA (jp-minerals.org/vesta/en/download.html)

References

For further information about the data processing in PETS2.0, see:

- L. Palatinus et al. Specifics of the data processing of precession electron diffraction tomography data and their implementation in the program PETS2.0. Acta Cryst. B 75: 512-522 (2019).

-PETS2.0 manual pets.fzu.cz/download/

For further information about the dynamical refinement and associated parameters:

- L. Palatinus et al. Structure refinement using precession electron diffraction tomography and dynamical diffraction: theory and implementation. Acta Cryst. A71: 235–244 (2015).

- L. Palatinus et al. Structure refinement using precession electron diffraction tomography and dynamical diffraction: tests on experimental data. Acta Crystallogr B71: 740–751 (2015)

Instructions

PART 1 - Data reduction using PETS 2.0

1. Data import

File B18H22.pts is the input for PETS. The file was generated automatically during data collection.

Open file B18H22.pts in any plain-text editor (Notepad, Notepad++, Vim, ... NOT MS Word). Lines beginning with a hashtag are ignored by PETS. The relevant lines are the following:

```
lambda 0.0251
Aperpixel 0.001850
dstarmax 1.2
phi 1.00
omega 0.5
noiseparameters 3.5 38
reflectionsize 20
I/sigma 5
bin 2
imagelist
dp-100\001.tif -20.00
                       0.00
dp-100\002.tif -19.00
                        0.00
. . .
endimagelist
```

Lambda: relativistic wavelength in Å Aperpixel: size of one pixel in Å⁻¹ phi: precession angle omega: orientation of the tilt axis (alpha) bin: binning factor by which the image resolution is reduced for image processing reflectionsize: the reflection diameter. Other parameters, especially concerning the detector specific noiseparameters, are explained in more detail in the manual available at http://pets.fzu.cz/ Open file Calibrated_optical_distortions.txt and copy the whole text into the B18H22.pts anywhere in the text with the exception of the imagelist. The resulting B18H22.pts file may look like the following:

lambda 0.0251 Aperpixel 0.001850 dstarmax 1.2 phi 1.00 omega 0.5 noiseparameters 3.5 38 reflectionsize 20 I/sigma 5 bin 2 #If yes, center of distortions is assumed to be a constant offset from the center of the diffraction pattern. Otherwise constant coordinates of the distortion centers are assumed. distcenterasoffset no #Distortion parameters.
 .ortions
 .040.0000
 990.0000
 0.0000
 0.00000

 0.000000
 0.000000
 1.100000
 0.000000
 0.000000

 0.000000
 0.000000
 0.000000
 0.000000
 0.000000

 0.000000
 0.000000
 0.000000
 0.000000
 0.000000
 distortions 1040.0000
 0.000000
 2.500000

 0.000000
 0.000000

 0.000000
 0.000000

 0.000000
 0.000000

 0.000000
 0.000000
 0.000000 0.000000 0.00000 0.000000 10.000000 0.500000 0.000000 0.000000 0.00000 0.00000 0.00000 enddistortions #Refinement keys for the distortion parameters. distortionskeys 0 Enddistortionskeys imagelist dp-100\001.tif -20.00 0.00 dp-100\002.tif -19.00 0.00 Endimagelist

The first two numbers in the first row of "distortions" are the coordinates of the optical axis of the microscope on the detector. The pincushion distortion represents the third number in the second row and the spiral distortion is the third number in the fifth row. Amplitude and phase of the elliptical distortion are determined in the rows four and seven by the first and fourth number, respectively. All values are in pixels. All distortion keys are zeroes, which means that the distortions will not be refined during the unit cell refinement. Keep your text editor opened, we will update some parameters later.

2. Run PETS 2.0

Start PETS by double clicking the executable Go to File->Open in the upper menu Open the file B18H22.pts in the folder

In the Main panel (Image data tab) you can view individual diffraction images (frames) using the arrows at the bottom of the panel. The contrast of the images can be adjusted by changing

the value of "Display cut off" at the bottom of the starting window showing the "Image data" panel. At e.g. 15000 only the primary beam and the strongest reflections are visible, at 10 we can clearly see the detector background noise. Values between 50 and 500 are useful to look at the relevant features, i.e. the reflections.

The arrows to the left of "Frame" (at the bottom controls below the diffraction pattern) and the value next to it control the image that is shown. The options "raw" and "processed" show you the raw diffraction pattern as recorded and the processed frame (some background correction etc.), respectively. "Simulated" frames will be available after frame orientation procedure.

3. Parameters



Data reduction steps are run by clicking the respective action button. The action options can be modified in the panel, which opens after clicking the small arrow next to the action button. The "Parameters" action button (first in the list of action buttons) is gray because there is no associated action. Note that you can adapt the size of the different window components (e.g. make the plot window more narrow or wider) and decouple the action bar or the console from the main window.

Click on the arrow (v) next to "Parameters" to check the following options:

In our input file .pts, we defined "reflectionsize 20" and "bin 2", hence the binned reflection diameter is 10 pixel. The default parameters used for min d* of 0.05 Å⁻¹ is OK. Here it means that any reflection/peak that is less than 0.05*Aperpixel pixel away from the center is ignored. "Precession" geometry was set in the input file because the data set was measured using precession.

4. Peak search

Click on the arrow (v) next to "Peak search", keep all setting and run "Peak search" by clicking on the action button.

The console output reports that 4178 peaks were written to the raw peak list (saved as B18H22.rpl). During peak search the determined pixel coordinates of the center of each diffraction pattern is determined and saved as the file B18H22.cenloc. Click on "Graph" (next to "Image data" and "Section images").

Red line: deviation of the x coordinate of the center from its initial value as a function of the frame number. Green line: deviation of the y coordinate. This should be checked because in the peak search options we used "Always accept center", which might fail. You can see a sawtooth shape of the lines, which is a result of a crystal shift within the convergent beam caused by the goniometer tilt, which is corrected every fourth frame. On top of that, after especially frames 28, 52 and 56 there are visible jumps of the center. These jumps are real and they are a result of

shifting the beam to the different (fresh) crystals spots as it is marked at the top of the figure. Last eight frames are taken from the same spot as the first 28 frames (marked as spot 1B). At the current stage, we assume that all pattern centers were correctly determined.



5. Rotation axis (omega)

The next step after peak search is the refinement of the angle omega between the horizontal axis (x axis of the diffraction image) and the projection of the rotation axis. This angle depends on the experimental setup, and needs to be refined for each data set. Click the action button "Rotation axis" to start the procedure.

In the plot window "Cylindric projection", the difference space of the extracted peak positions is represented as a cylindric projection. In the ideal case, it consists of well-defined spots on sinusoidal curves and no elongated features. Here, the starting angle (defined in the input file) is 0.5° and it refines to 0.1562.

6. Peak analysis

Click the action button "Peak analysis".

PETS starts the analysis of the peak positions obtained in the peak hunting procedure. In the first step, distance distribution between peaks in the image plane is analyzed. A sorted plot of inter-peak distances (left screenshot) is displayed in the Processing tab (red curve) with its derivative (green):



For good quality data sets, the red curve contains distinct steps, and, consequently, the green curve has distinct peaks.

Peak analysis(continue) Click "Peak analysis (continue)"

Difference peaks are saved to B18H22.diff

In the next step, an autoconvolution of the diffraction pattern (difference space) is analyzed, and the groups of peaks in the autoconvolution (clusters) are replaced by the cluster centers. Again, a distance distribution is displayed (right screenshot), and clear-cut steps on the curve indicate a well-defined lattice.

Click "Peak analysis (continue)"

Cluster map file is written to B18H22.clust

Now the first step is finished and PETS prepared the files needed for the determination of the orientation matrix and the unit cell parameters.

7. Find unit cell and orientation matrix

Click "Find unit cell and orientation matrix"

The main menu on the left is replaced by a new set of buttons to determine, refine and modify the orientation matrix, which in turn defines the unit cell parameters. On the right side, a new tab "3D panel" opens, where it displays a projection of the peaks determined from the peak search and peak analysis. Note that the view can be enlarged by scrolling the wheel of your mouse, and rotated by holding the left mouse button. By default, the "Data used for indexing" is "xyz", i.e. the points in B18H22.xyz are used for the plot and for the orientation matrix determination.

Click "Find possible cells automatically"

Data used for indexing: vyz v	
Cell 1 indexed/all: 1381 / 1499 indexed: 92.13 % volume: 1390.57 Bravais class: tP	
New cell Delete cell	
Undo Redo	
a b c α β γ cell: 10.7474 10.7790 12.0037 89.929 89.824 89.877 s.u.: ⊙ Orientation matrix	
-0.052996 0.076360 0.001051	
0.046448 0.033054 -0.065925	
-0.000739 -0.041050 -0.030922	
 Find possible cells automatically 	
(a) from difference space (b) from triplets Maximal d* difference for indexing (rec. b): (0.018)	
Angular tolerance for rect direction (deg):	
Maximai Volume (A^s):	
Symmetry search: tolerance for lengths (%): 3	
tolerance for angles (deg): 0.6	
a b c α β γ V Bravais	ind/all

Open the details by clicking on the arrow next to the action button. You see the parameters used in the procedure as well as the possible unit cells found by PETS (only one is found with the default settings). The first one in the list (see screenshot) is automatically used. The cell angles are very close to 90 degrees, a and b are almost equal. Therefore, PETS determined that the Bravais class is tetragonal primitive (tP). Based on the initial orientation matrix, reciprocal space sections can be generated to investigate visually the symmetry in reciprocal space. As a part of the automatic indexing procedure the lattice parameters and orientation matrix are refined. We are happy with the result and thus no additional cell refinement is needed at this moment.

Leave the indexing panel by clicking "Finish"

8. Create reciprocal-space sections

In the parameters (arrow next to action button "Reciprocal-space section"), you can change the size of pixels in the reconstruction images and the thickness of the slab that is projected into the section in $Å^{-1}$. Too small numbers may result in "holey" reconstructions. We will proceed with the default values.

Click "Reciprocal-space sections"; OK (use the standard set)

PETS reconstructs sections through reciprocal space at planes hk0, hk1, hk2, h0l, h1l, h2l, 0kl, 1kl, 2kl from the diffraction patterns. These sections are used as an aid to determine the space group. The generation of the sections takes about 2 minutes. Change to the tab "Section images" Choose "hk1.tif" Set the Display cut off to 30 The hk1 section, though very incomplete, clearly violates tetragonal symmetry. The point group of this section seems to be 2mm. Note that the white circle with a cross in the middle of the image indicates the center of the section, in this case the expected position of the 001 reflection. In many cases, the reliable determination of the Laue class of the tetragonal or cubic crystal systems requires checking sections that are not included in the standard set, e.g. sections perpendicular to <110>.

The sections hk0, h0l and 0kl may give indications on the presence of glide planes and/or 2_1 screw axes. For your convenience you may turn on the grid, which will show the positions of the expected diffraction maxima



by clicking on "gr" button next to the Display cut off window. The Okl section suggests the presence of an n-glide plane as reflections with k+l = 2n+1 are mostly absent or very weak. The hOl section suggests the presence of a-glide plane (h odd are either absent or very weak) perpendicular to the b-axis. Only a small fraction of the hk0 section was measured. However, this is sufficient to recognize that reflections with h odd are absent, suggesting a-glide plane perpendicular to the c-axis. From this quick analysis we derive that the space group Pnaa most likely describes the symmetry of the crystal structure. The standard setting of this space group (#57) is Pccn and we should therefore adapt the orientation matrix.



9. Change orientation matrix

Click "Find unit cell and orientation matrix" Click on "..." next to the Bravais Class Change the Crystal system to "o" (orthorhombic); OK The Bravais class is now "oP" Open the parameters of "Modify cell" (click on the arrow)

Data used for indexing:	хуz	~		
Cell 1	indexed/all: 1381 / 1499 indexed: 92.13 % volume: 1390.57 Bravais class: oP			
New cell Delete cell Undo Redo a b cell: 10.7474 10.7790 s.u.: - - Orientation matrix - 0.052996 0.076360 0.046448 0.033054 - -	c α β 12.0037 89.929 89.824 89.8	1 777 -	$\begin{tabular}{ c c c c } \hline & & & & & & & & & & & & & & & & & & $	•

Enter the transformation matrix (right screenshot) 0 1 0; 0 0 1; 1 0 0; Click on "Transform by matrix"

The new lattice parameters are 10.7790 12.0037 10.7474 89.824 89.877 89.929. Important! Now the unit cell will be refined and for the best result it is neccesary to switch the data from .xyz to .cor, which contain the searched peaks aligned using the diffraction centers. Using .cor file, the refinement will be done against the data, which were not matematically processed in the step "Peak analysis".

Change "Data using for indexing" to .cor.

Open the parameters of "Refine cell" (click on the arrow).

Select "refine cell and distortions", check "refine cell". Leave all other settings default.

^				Refine c	ell		
⊖ refine UE ⊖ refine ce ● refine ce	l + cell Il from d Il and distortion	s					
Maximal d*	Maximal d* difference for indexing[Inv(A)]: 0.018						
symmetric	у		tri	clinic		~	
O fixed cel	parameters						
10 10 1	0 90 90 90						
fixed	up to a scale fa	ctor .					
✓ refine ce	II refine di	stortions					
Distort	ons						
Center:	x	1040					
	У	990					
rad	lin	quad	c	ub	phase		
0	0	0	0 1.1				
1	0	0	0	0	0		
2	0.5	0	0		55		
sg term:	s 🗖 0	0	0				
tan	lin	quad	c	ub	phase		
0	0	0	2.5				
1	0	0	0		0		
2	0.5	0	0	0	10		
sg term:	5 🗌 0	0	0				
🗌 tie el	liptical distortio	n					

The distortion parameters, which were copied at the beginning to the B18H22.pts can be seen by clicking on the "(v) Distortions" button. Click on "Refine cell"

	а	b	с	α	β	γ
cell:	10.7826	12.0118	10.7441	89.831	89.885	89.869
s.u.:	0.0009	0.0028	0.0016	0.023	0.010	0.016

The resulting unit cell angles are all by more than 0.1° off the expected 90°. This deviation is caused solely by a deterioration of the crystal, which results in a gradual change of the lattice parameters with the accumulaed electron dose as we will see in the part dedicated to frame-by-frame distortions refinement. To mitigate this problem, we will use the symmetry restrictions of the unit cell during the refinement.

Select the "symmetry" orthorhombic and click on "Refine cell".

^		Refine cell
00	refine UB + cell refine cell from d refine cell and distortions	
M	aximal d* difference for indexing[Inv(A)]:	0.018
۲	symmetry	orthorhombic v
0	fixed cell parameters	
	10 10 10 90 90 90	
	fixed up to a scale factor	
	refine cell refine distortions	

The resulting lattice parameters are 10.7826 12.0118 10.7441 90.000 90.000 90.000 Click on "Finish"

The orientation matrix is written to B18H22.celllist. Note that the reciprocal space sections were not changed and thus do not correspond to the new orientation matrix. You should generate new sections.

Click "Reciprocal-space sections"; OK (use the standard set)

10. Frame orientation optimization

Open the options of "Optimize geometry and integration parameters".

Uncheck the "refinement using rocking curve profile" and check "refinement using frame simulation".

Optimize geometry and integration parameters						
refinement using rocking curve profile (camel plot)						
🗹 RC width (rec. Å)	0.001					
🗹 apparent mosaicity (deg.)	0.05					
precession angle (deg.)	1					
calibration constant (rec. Å, binne	d) 0.0037					
minimum l/sigma	10					
 ✓ refinement using frame simulation ✓ orientation angles ✓ center of the diffraction patterns ■ C width ✓ mosaicity → distortion for simulation use: 	☐ Transfer to global ④ uniform intensities ◯ integrated intensity ttings					

Click on "Optimize geometry and integration parameters".

The procedure will use default uniformly decreasing intensities with increasing resolution. PETS will find optimal values of the three correction angles alpha, beta and omega.



You can check the process in "Image data" card by switching the radiobutton to "simulated". Then you can choose if you want to see the data and the simulation in parallel or if you fant to see the combination of the two images (purple – data, green – simulation). Or if you want to see just the simulation.



^

Check the "distortion".

^ Optimize ge	Optimize geometry and integration parameters				
refinement using rocking curve profile (camel plot)					
RC width (rec. Å)	0.001				
apparent mosaicity (deg.)	0.05				
precession angle (deg.)	1				
calibration constant (rec. Å, binned)	0.0037				
minimum l/sigma	10				
 ✓ refinement using frame simulation ✓ orientation angles ✓ center of the diffraction patterns ○ C vidth ✓ mosaicity ✓ distortion for simulation use: 	Transfer to global Guniform intensities integrated intensity				

This will allow PETS to refine individual frame-by-frame magnification and amplitude and phase of the elliptical distortion.

Click on "Optimize geometry and integration parameters".

The correction angles are now much smaller.

Choose "Frame-by-frame" distortions in the Graph tab.

There you will see the correction of the apparent magnification (red curve), amplitude and phase of the elliptical distortion (green and blue curve, respectively). The elliptical distortion is constant, which is normal for the microscope setting, which was used during this data collection. In terms of the lattice parameters stability this indicates that the lattice parameter ratios and angles remain constant. This should be the case also for the apparent magnification (red curve), if the unit cell volume woulb be constant. But there is an obvious trend observable in the graph which is strongly correlated with the accumulated doses in the particular measured spots. The apparent magnification of the reciprocal space is decreasing, which means that the actual volume of the unit cell increases as a result of the accumulated beam damage.



11. Integrate intensities

Click on "Process frames for integration" and uncheck "shift integration box to maximum intensity".

All the significant distortions in the data were corrected so the integration box shift is no more necessary during the integration. In the "Image data" tab, diffraction images with circles (integration boxes) drawn at positions predicted from the orientation matrix are displayed. Green circle: observed reflection with I>3 σ (I). Brown circle: reflection with I<3 σ (I). Diamond: reflection expected either on upcoming or on previous frame(s). In the parameters the "max d* for integration" was set to 1.2 Å⁻¹. The main reason for this is that a reduced resolution used for the refinement will reduce the time needed for dynamical refinement cycles later. The curious user is encouraged to set the "max d* for integration" limit to 1.4 in the options of "Parameters".

In the "Graph" tab, choose the "Rocking curve"

From these plots the RC width and the apparent mosaicity can be derived.

Open the options of "Optimize geometry and integration parameters"

Note the default options: RC (rocking curve) width = -0.001 Å⁻¹, apparent mosaicity = 0.05°. Check the "refinement using rocking curve profile" and uncheck "refinement using frame simulation".

•	Optimize geometry and integration parameters						
	refinement using rocking curve profile (camel plot)						
	🗹 RC width (rec. Å)		0.001				
	🗹 apparent mosaicity	(deg.)	0.05	1			
	precession angle (d	eg.)	1	1			
	calibration constan	t (rec. Â, binned)	0.0037				
	minimum l/sigma		10				
	refinement using frame orientation angles center of the diffrac RC width mosaicity distortion for simulation use:	: simulation tion patterns Advanced settings	☐ Transfer to global	is ty			

Click on "Optimize geometry and integration parameters".

Note how the parameters change once the parameters are refined. The predicted camel plot (red curve) fits now better to the observed profiles (blue curves).

Left screenshot: before optimized geometry. Right screenshot: with optimized geometry The determined RC with is 0.0006 Å⁻¹ and the apparent mosaicity 0.032°.



12. Intensity integration

Open the parameters for "Finalize integration". Check that both "kinematical" and "dynamical" integrations are checked. Activate "Frame scaling".

Set the Laue class to "mmm" (screenshot).



Click on "Finalize integration"

In the console, the Rint is shown for different Laue classes. If we finalize integration with two different "intensity estimation" settings ("fit profile" and "integrate profile"), we end up with the following statistics for the Laue class mmm:

Intensity estimation	Rint(obs)	Rint(all)	Nobs/Nall
Integrate profile	20.39	21.39	830/1251
Fit profile	12.50	13.91	729/1242

Default value for intensity estimation is "integrate profile", which is the more robust technique and works well even for poorer data. However, for good data the "fit profile" option gives better results.

Select the radio button "fit profile" for the "intensity estimation" Click on "Finalize integration"

Two files were saved in the main folder: B18H22.cif_pets and B18H22_dyn.cif_pets. Look at the "Frame scales" of the Graph tab



13. Unit cell refinement after the removal of the effect of the beam damage

Click on "Find unit cell and orientation matrix". Open "Refine cell" options.

Set "Maximal d* difference for indexing to 0.01 A⁻¹ and "symmetry" to triclinic. Click on "Refine cell".

The resulting lattice parameters are 10.785(1) 12.012(2) 10.744(1) 89.99(2) 89.98(1) 89.98(1). The lattice constants obtained from XRPD data measured at the same temperature were 10.7566(3) 11.9796(4) 10.7142(3) 90 90 90. The a, b, and c lattice parameters from electron diffraction are 0.26, 0.27, and 0.28% longer than those from the XRPD. The source of this discrepancy is twofold. The first is the unit cell expansion, which is dependent on the accumulated dose. From the image below, you can see that the first four frames taken on the fresh spots have on average the apparent magnification by 1.0 pixels larger than the magnification represented by the refined unit cell. Therefore, the unit cell parameters on these frames are 0.19% shorter than the parameters of the unit cell obtained in step 9. The resting ~0.08% difference comes from the unevitable fluctuation of the magnification, which according to our experience is about 0.2% for the excitation of the diffraction lens close to the eucentric focus. The fluctuation is not reflected in the e.s.d.s. In conclusion, using calibrated microscope

distortions and beam damage effects suppresion it is possible to arrive to the lattice parameters very similar to those obtianed by the XRPD.



Click on "Finish". You can close PETS now. Answer "Yes" to the question about saving the project

The current state of the project is saved in the file B18H22.pts2. Opening this file with PETS next time allows you to continue from the current state of processing. If you open B18H22.pts, you must start from scratch again, because B18H22.pts is not overwritten by PETS and still stores the initial settings.

PART 2: Structure solution and kinematical refinement in Jana2020

1. Create new structure

Create a new folder with the name "Jana_kinematical"

Important! The data-processing procedure is almost never perfectly reproducible. Small differences in the indexing and cell refinement procedure may result in small differences of integrated intensities. If you want to be sure that you can reproduce the following part of the tutorial, do not use the file B18H22.cif_pets that you just created, but use the file B18H22.cif_pets from the folder PETS_output_CIF. Using your own cif_pets file is also possible, but your results may slightly differ from the results described in this tutorial. Copy the file B18H22.cif_pets to the folder "Jana_kinematical" Start Jana2020 In the Main menu bar, use "Structure → New" and open new structure "B18H22" in directory

2. Import Wizard

Jana_kinematical

[On the screen: Specify type of the file to be imported] Select "Single crystal: Known diffractometer formats"; NEXT [On the screen: Data reduction file from] Select "Pets electron diffractometer" (Shortcut ALT + P) File name has automatically changed to B18H22.cif_pets. NEXT

Set temperature to 100; NEXT

Note the wavelength of 200 kV electrons, which is 0.0251 Å. Leave all settings unchanged (no twinning, no modulation); NEXT The program reads 3690 reflections from the file OK

For absorption correction select "None or done before importing"; NEXT FINISH; OK to accept made changes, i.e. accept the new data set

3. Symmetry Wizard

NEXT to close the information window and start the symmetry wizard.

Symmetry Wizard can be started separately by expanding "Reflection file" in the Command tree and selecting "Make space group test". The default settings are prepared for x-ray diffraction. For electron diffraction we have to increase the tolerances.

Set "Maximal ave(I/sig(I)) for systematic extinctions induced by other symmetry" to 25 (default is 10)

Thus, less space groups will be excluded even if there are many intense reflections in the data set that are forbidden by the reflection conditions of the respective space group.

NEXT;OK

Select Orthorhombic Laue symmetry (mmm) NEXT

The statistics clearly exclude a centered cell. The choice of a primitive cell is also in agreement with our visual inspection of the reciprocal space sections.

NEXT

The reciprocal space sections suggest space group Pccn. The intensity statistics indicate that 41 observed reflections with an average $I/\sigma(I)$ of 11.487 contradict the reflection conditions of Pccn. For electron diffraction experiments, this may be acceptable and a careful analysis of the reciprocal space sections often results in a correct determination of the space group. In our case, we have good indications for the n glide plane perpendicular to c^* and the c glide plane perpendicular to a^* . Pccn is the most probable space group complying with these observations. Choose Pccn; NEXT

Accept the space group in the standard setting; FINISH Symmetry is saved in file B18H22.m50.

4. Creating refinement reflection file

In this step, the program creates the file B18H22.m90 containing the data set merged by symmetry and with discarded forbidden reflections. This file will be used for refinement. Proceed with the wizard, leaving all settings default.

The summary after averaging indicates an internal R-factor of 12.34 for 699 observed reflections.

OK; FINISH

5. Structure Solution Wizard

OK to close the information window and start the solution wizard.

Structure solution wizard can be executed separately using the quick button, or from the Command tree through "Structure solution->Run Superflip".

30°	
\odot	Phases
, 🕄	Edit structure parameters
47	Structure solution
L P	Run Superflip
5	🗄 Fourier synthesis

In "Formula" textbox type: B9 H11 In "Formula units" type 8 Check "Use a specific random seed" and enter the value 2. Switch iteration scheme to AAR. "Run Superflip".

Commands for Superflir								
Commanus for Supering								
Basic commands Ad	vanced commands							
Formula:	B9 H11		Phase:					
Formula units:	8 Calculate density	Sum formula from st	ructure model					
Actual space group:	Actual space group: Pccn Change the space group							
🗌 use in le Bail dec	omposition structure information for alrea	dy identified phases						
allow manual edit	ting of the command file before start							
use previously pr	epared input file for Superflip							
use old solution a	and reinterpreted							
Repeat Superflip	: Until the convergence detected		Biso: 0					
🖂 Repeat Superflip	: Number of runs => 10		Maxcycles: 2000					
🗌 Use local normaliz	zation							
Use a specific ran	ndom seed => 2							
Define explicitly of	delta value => 0.9							
Iteration scheme:	⊖ CF F	or peak search use:	O EDMA - fixed composition					
	OLDE		○ EDMA - fixed number of atoms => 0					
	AAR		\bigcirc EDMA - peak interpretation by Jana2020					
Starting model:	Random phases		Peaks from Jana2020					
	O Patterson superposition map		O Peaks from Jana2020 but first run Fourier					
	Run Superflip Open the	listing Draw str	ucture Draw 3d map					
	Use tutorial files Accept last	solution Qu	it					

The specific random seed is used just to guarantee reproducibility of the results. This also requires that the imported data set (B18H22.cif_pets) is identical to the one used for this tutorial.

Window of Superflip appears with iteration running. It finishes after 10 runs. Check "Use tutorial files"

This option ensures the reproducibility of the structure solution step. Jana2020 will replace the job. Files by job.tut_files, which were generated in the same way as it is described in this manual. However, due small randon effects, which occur during the structure solution, small differences were frequently found in between the solutions. Therefore, we have decided to provide the solution files to ensure the reproducibility. You can of course use your own solutions, but the agreement between your results and the results in this manual may slightly differ.

Information about the match between the found density symmetry and the suggested space group is shown after the end of the search.

)	INFORMATION	
	R-value after noise supressio HM symbol: Pccn	on: 31.332
1	Symbol R-va	lue
	n 0.82	2
Ĩ	с 8.04	1
	c 8.27	7
1	Formula: B9 H0	
-	-	1
	OK	
f	-	

OK; for reproducibility reasons click on "Replace the result with tutorial files"; OK. You can also go back to your original result using the same button.

Run Superfli	Open the listing	Draw structure	Draw 3d map
Replace the result with tutorial fil	Accept last solution	Quit	

Press "Open the listing"

Superflip performed ten solution attempts, and picked the best one. Its characteristics are displayed near the bottom of the listing:

At the very bottom the most reliable solution is described after the line "Properties of the saved density":

Run	Rvalue	Peaks	Symm.	Der.SG
1	30.54	3.00	8.67	Pccn

The space group proposed by Superflip is the same as our initial guess. By going through the complete log file we can see thatfor all solutions out of the 10 attempts the space group was Pccn, further supporting our initial analysis of the space group symmetry. We can now close the log file. Note that it is saved B18H22.sflog, so you can also open it with a text editor. Close the listing

Click on "Draw 3d map"

VESTA opens with a density shown as isosurfaces (in the default settings the shown isosurface level corresponds to the 2 σ level based on the statistics of the 3D electrostatic potential map). The map is shown together with an interpretation of it as a model of boron atoms (green spheres). Note that if you used a different B18H22.cif_pets file, most likely your solution will look different. The shape of the blobs will actually not differ much, but an origin shift (most likely 1/2a, 1/2b, 1/2c or a combination of these) will make a direct comparison more challenging. The solution puts the molecule in the middle of the unit cell. You can limit the 3D map by setting the boundaries on the card "style" to the limits shown in the figure. Tools (from the menu bar) ->Boundary (CTRL+Shift+B); change the range of x, y, z as in the screenshot; APPLY; OK



Set viewing along b (click on "b" in the upper left corner).

6. Hydrogen atoms in the structure solution

The last view in VESTA shows no signs of any of the 11 hydrogen atoms. However, changing the isosurface level to 0.4 (approximately 0.75σ) reveals new local maxima in the density map. In VESTA:

Tools->Properties->Isosurfaces (CTRL+Shift+P); click on "New"; change the isosurface level to 0.5; mode "Positive"; change the color to white; OK



Most of the new local maxima visible at the lower isosurface level are close to the identified boron positions. Actually, 10 of them correspond to the hydrogen positions (11 total) we are looking for. Close VESTA (do not save the file)

 Run Superflip
 Open the listing
 Draw structure
 Draw 3d map

 Use tutorial files
 Accept last solution
 Quit

In Jana: keep checked the "Use tutorialfiles" and click on "Accept last solution".

The program creates file B18H22.m40 with 9 B atoms.

The state of the Jana files after clicking on "Accept last solution" can be found in the folder "JANA/O-solution/" You can work with the files in this folder if you prefer. In that folder, the B18H22.sflog is the log file of the structure solution and the M40 file contains the resulting structural model (9 B atoms, list of 72 local maxima). Opening this job with Jana you can open the electrostatic potential map with Run Contour->Use old maps->draw maps as calculated.

7. Kinematical refinement and identification of the hydrogen atoms

[On the screen: basic window of Jana]

Double click on "Refinement"-> "Run refinement"

The refinement converges after 7 cycles to R(obs) 24.30. At the end of the refinement Jana notifies you that there are serious warnings, which are related to (physically meaningless) negative displacement parameter of the B2 atom. We now want to look at the difference Fourier map to identify the hydrogen atoms.

RFactors	overview
R factors GOF(obs R(obs)= Last wR Maximun	s: [1127=703+424/37], Damping factor: 1.0000 s)= 14.14 GOF(al)= 11.11 24.30 wR2(obs)= 51.74 R(al)= 28.12 wR2(al)= 51.96 (al): 59.13 56.74 52.01 51.95 51.96 51.96 51.96 n change/s.u. : -0.0284 for y[85]
	Question
	Serious warnings in the listing! Regular end of REFINE program
	Yes No

Double click on "Run Contour"

"Calculate a new Fourier map" -> "draw maps as calculated" In the card "Basic", select "Map type:" "F(obs)-F(calc) – difference Fourier"

	□- Calculate a new Fourier map draw maps as calculated		
F	ourier commands		
	Basic Scope Peaks		
	Map type: F(obs)-F(calc) - difference Fourier		

In the page "Scope", select "Whole cell"; OK

A contour map is shown, which covers the coordinate range $0 \le x \le 1$, $0 \le y \le 1$, z = 0However, the identification of relevant peaks in the Difference Fourier map is easier with a program supporting a 3D view.

In the menu bar, use "Run->Vesta"



Alternatively, VESTA can be started from the lower toolbar

Objects->Boundary (CTRL+Shift+B); change the range of x, y, z as in the screenshot; APPLY; OK Ranges of fractional coordinates

x(min) =	0.1	x(max) =	0.9
y(min) =	0.2	y(max) =	0.8
z(min) =	0	z(max) =	1

Objects->Properties->Isosurfaces (CTRL+Shift+P), Isosurfaces Click on the first entry in the list (1, 0.220, Positive and negative) Enter in the input field of "Isosurface level": 0.220*1.25 (This will set the isosurface level to 2.5σ)

Isosuraces				
F(min) = -0.415070; F(max) = 0.535830;				
Render from front to back				
Positive and negative V Opacity 1 (0~255): 127				
Isosurface level: 0.274686 Opacity 2 (0~255): 255 🔹				
Color: 255 🐳 255 🐳 0 💌				
No. level mode color New				
1 0.274686 Positive and negative Delete				

All the 11 hydrogen atoms are clearly visible.

Close Vesta "Quit" the Contour program in Jana

Right click on the quick button "Run Fourier synthesis"

Check that the Map type is "F(obs)-F(calc) – difference Fourier"

In the section "Peaks": Maximum number of peaks, select "Define them explicitly", Positive: 30

Basic	Scope	Peaks							
					🗹 A d	tivate	peak sea	arch procedure	1
						Maxin	num nun	nber of peaks:	
(🔾 Use de	efault va	lues						
(Define	them e	xplicitly	Po	ositive:	<u>β0</u>	•	Negative:	5

OK; "Yes+Start" to start the Fourier calculation.

YES to start the procedure of including of new atoms. Set the minimal and maximal distance as on the screenshot.

Skip peaks being too close to existing atoms		ting atoms	Show distances up to:	1.5	Angs.
Minimal distance:	1	Angs.			

The peak list is shown with the option to "Include the peak at the specified position" as new atom.

Include Max1 as H1b3b8 Include Max2 as H1b3 Include Max4 as H1b2 Include Max7 as H1b7 Include Max8 as H1b4 Include Max10 as H1b3b7 Include Max12 as H1b8 Include Max13 as H1b5 Include Max15 as H1b1b9 Include Max24 as H1b6

Include Max33 as H1b1

We have now 9 B atoms and 11 H atoms, i.e. all expected atoms were identified. We can now start the final refinement.

Note: to verify the proper number of atoms, start the quick button "Edit basic paraneters", go to the page "Composition" and click "Sum formula from structure model". It will show B9 H11 for eight formula units.

Finish; YES.

Close listing of Fourier calculation.

Double click on "Run Contour"->"Use old Fourier map"->"draw map as calculated" and run Vesta from the lower tool bar. Set isosurface level to 2.5σ and boundaries as in the screenshot. Check the positions of the introduced hydrogen atoms.

Close Vesta, close Contour.

Double click on "Refinement -> Run refinement" to start the refinement.

RFactors ov	rerview
R factors : GOF(obs)= R(obs)= Last wR2(a Maximum c	[1127=703+424/81], Damping factor: 1.0000 • 8.70 GOF(al)= 6.78 15.84 wR2(obs)= 32.77 R(al)= 19.96 wR2(al)= 33.07 al): 37.91 33.34 33.16 33.10 33.09 33.08 33.07 thange/s.u.: -0.0455 for z[H1b8]
	Question
	Regular end of REFINE program
	Open the listing?
	Yes No
	It remains 6 sec - press this button to continue

The refinement converges after 7 cycles.

Refinement -> View Refinement listing -> Go to (ALT+g) -> BOTTOM and scroll up. There are no "Serious warnings" so the displacement parameters of all atoms are positive. Close listing.

	 Difficultionshops and fear 1 = fear and 1 = 0
📕 Jana20)20 - [D:\Workshops Jana\2021_2\Borane\Jana_PB\B18H22
Structure	Transformations Settings Tools Help
Main	
	■ ♣

Click on "JanaDraw" and in the lower toolbar on "options" — and set the "Bonds between

atoms...expanded by:" to 20% and then click on "Build molecule"

The displacement parameters can be visualized in JanaDraw by clicking the icon $^{\textcircled{}}$ (Draw ellipsoids).

The structure looks alright and displacement parameters seem reasonable. However, the ADPs cover rather large range of values. As a rule of thumb, the same atoms in a similar environment should have a similar (equivalent) isotropic displacement parameter.



Close JanaDraw.

Double click on "Distance calculation" and "Distance commands". Uncheck "d(max) derived from atomic radii and set the d(max) value to 1.7.



OK;Yes+start;Yes. Check the B-H distances in the listing.

H1b1b9, H1b3b7 and H1b3b8 form 3-center 2-electron bonds and their bond lengths range from 1.21 to 1.38 Å. For the other hydrogens the bond lengths range from 1.16 to 1.46 Å. The uncertainties are between 0.04 and 0.05 Å. We can now move on to the dynamical refinement.

The Jana files at the end of the kinematical refinement can be found in the folder "JANA/1kinematical" for comparison. The following part of the tutorial is based on the files from this refinement.

PART 3: Dynamical refinement in Jana2020

1. Preparation of reflection file for the dynamical refinement

During the integration procedure in PETS, an input file for dynamical refinement was created. The integration for the dynamical refinement considers each diffraction image separately. The result is a list with the hkl indices, the measured intensity with the respective sigma and the number of the frame in which the intensity is present. PETS creates the file B18H22_dyn.cif_pets with all this information, which will be used in Jana2020.

Create a new folder "Jana_dynamical"

If you want to be sure that you can exactly reproduce the following part of the tutorial, use the file B18H22_dyn.cif_pets from the folder PETS_output_CIF provided with the example files. Copy the PETS output file "B18H22_dyn.cif_pets" into the folder Jana_dynamical [On the screen: basic window of Jana after finishing the kinematical refinement] "Structure→Save as"

Locate the folder "Jana_dynamical" and save a new structure B18H22_dyn Continue with the new structure.

"Reflection file \rightarrow Import/modify reflection file \rightarrow Reimport"

Select the B18H22_dyn.cif_pets as the new reflection file.

Make the reflection file for dynamical refinement

Check "Make the reflection file for dynamical refinemant; NEXT.

Set temperature to 100; NEXT; NEXT.

Jana reads 16707 reflections.

OK; NEXT; FINISH; OK.

Process the reflection by clicking on NEXT.

Notice that 4956 reflections are observed out of 16707 reflections read. 79 observed reflections were rejected because they violate the systematic extinctions.

OK;OK;FINISH.

"Edit structure parameters" → "Edit parameters for electron diffraction"

Set "Maximal diffraction vector g(max)" to 1.25

Set "Number of integration steps" to 200

Set "Limit on RSg" to 0.8 (this increases the number of reflections used for the refinement) Select "For Fourier rescale to Fcalc"

Set "Number of threads" to the number of cores of your processor.

Activate "For Fourier rescale to Fcalc".

Orientation matrix: U11 0.07634 U12 0.00119 U13 -0.05284 U21 0.03309 U22 0.06585 U23 0.04627 U33 -0.04094 U32 -0.05091 U33 -0.06107				
Maximal diffraction vector g(max): 1.25 Number of integration steps: 200				
Dyngo commands:				
Select zones for refinement Define zones of equal thicknesses				
Run optimizations except of scale, optimize also: Thickness Show thickness plots Orientation Orientation Orientation Orientation Orientation				
Zone#: 1 - uvw: (-0.677,-0.351,-1.000) R(al)= Select zones for editing				
Precession angle 1 Alpha -19.752 Beta -0.114				
Esc OK				

2. Determination of the thickness

In "Edit parameters for electron diffraction" check "Thickness" and then click on "Run optimalizations"

"Show thickness plots"

For Zone#1, the best R factor is observed for a thickness of 840 Å. But for Zone#3, ther are two minima about 500 and 1500 Å. Look at the rest of the frames. There are two frequent best thicknesses – one about 1000 and the other about 1500 Å. It will turn out during the refinement that the better choice is 1500 Å by 0.7 percent point of R(obs). If the thickness is larger than 1000 Å and the number of integration steps is about 100 you should consider to incease it to at least to 150. Thicker crystals have narrower diffraction maxima, which can be missed by a too low number of integration steps.



Close the plot.

Click on "Select zones for editing", "Select all", OK.

The EDThick box is yellow and locked because the parameter is not the same for all selected zones.

Click EDThick box; click Unlock and change the value to 1500. This changes EDThick for all zones to 1500 Å.

Click on "Select zones for editing"; "Refresh"; OK

By changing number of zone by "Zone#" textbox, we can see that all zones have EDThick 1500 Å Uncheck the checkbox "Thickness"

Click on "Run optimalizations"

This will optimize the scale factor for each frame based on a thickness of 1500 Å. The optimization runs for about 2 minutes.

Zone#: 1 • uvw: (-0).677,-0.351,-1.000)	R(all)=	Select zon	nes for editing
Precession angle	Alpha -19.752	Beta -0.114		
EDScale 1071.819	ED Thick 1500	EDPhi 0		EDTheta 0
Unlock	Esc	ОК		

OK; YES to save m42 file.

The current state of the Jana files (before clicking "Refine") is available in the tutorial folder: JANA/2-dynamic/B18H22_dyn_0-setup

3. Hydrogen localization

Click on "Edit atoms" icon 🖾 . Select all hydrogen atoms and delete them by clicking on

Action and selecting delete atoms.

Yes to all; OK; Yes for rewriting m40.

Run refinement. It converges after 5 cycles.

The resulting state of the Jana files is available in the tutorial folder as JANA/2dynamic/B18H22_1-hydrogens. You may use these files to save time, look at the results or simply for comparison purposes.

Double click on "Fourier synthesis" to unfold the options.

Double click on "Calculate Fourier maps"; No;No.

Click on "No"

Double click on "Run Contour"; double click on "Use old Fourier maps" and "Draw maps as calculated".

Click on Vesta icon ᆇ in the lower tool bar.

Objects->Boundary (CTRL+Shift+B); change the range of x, y, z as in the screenshot; APPLY; OK Ranges of fractional coordinates

x(min) =	0.1	x(max) =	0.9
y(min) =	0.2	y(max) =	0.8
z(min) =	0	z(max) =	1

Check difference maxima. The hydrogens are nicely visible even at 3σ level.

Close Vesta, close Contour.

Double click on "New">"New atom". Set limits for the maxima distance from the atoms as in the case of the kinematical refinement. All the first 11 maxima belong to the hydrogen atoms. Name the hydrogens in the same way as for the kinematical refinement (H1 + name of the boron atom(s) to which it is attached).

You can see, that dynamical refinement substantially improves the visibility of the difference maxima due to missing hydrogen atoms in the structure.

🖂 Skip peak	s being too close to existing atoms	Show	v distances u	ip to: 1.7	Angs.	
Minimal dista	ance: 1 Angs.					
	List of peaks	Equivalent co	ordinates of	the peak	Shortest dista	nce to
Max1	Max30	0.286467	0.153799	0.180884	as typed in	
Max2	Max31	0.713533	0.653799	0.319116	B7	1.226
Max3	Max32					
Max5	Max33					
Max6	Max35					
Max8	Max36					
Max10	Max37					
Max11	Max38					
Max12	Max39					
Max14	Max40					
Max18						
Max22						
Max24						
Max26						
Max28						
Max29						

Confirm the changes.

4. Refinement

Double click on "Refinement->Refinement commands", set the number of cycles to 20. OK; YES+start to start the refinement.

The refinement converges after 20 cycles to R(obs) 8.72% (3154 observed reflections) and wR2(all) 19.29% (8449 reflections). The total refinement time is about 2 hours (desktop computer, 24 physical cores). The resulting state of the Jana files is available in the tutorial folder as JANA/2-dynamic/B18H22_2-refine. You may use these files to save time, to look at the results or for comparison.

Close the output.

5. Optimization of orientation

Save a new structure "File \rightarrow save as \rightarrow B18H22_dyn_ori"; OK; Yes. "Edit structure parameters \rightarrow Edit parameters for electron diffraction". Enable "Except of scale optimize also: Orientation".

Run optimizations except of scale, optimize also:	□ Thickness Show thickness plots ☑ Orientation
Zone#: 1 • uvw: (-0.677,-0.351,-1.000)	R(all)= Select zones for editing
Precession angle 1 Alpha -19.752	Beta -0.114
EDScale 957.4786 D EDThick 1639.69	EDPhi 0 EDTheta 0

Click "Run optimalizations"

For each frame the orientation will be optimized.

Check the EDTheta parameter for each zone using arrows next to the zone number Only small number of frames have the orientation correction above 0.1°, which indicates that the frame orientation optimization done in PETS 2.0 produces precise frame orientations even for molecular crystal, which has only moderate resolution and deteriorates during the measurement.

Click on "OK" to leave the electron diffraction menu Click on "YES" to save the changes to the M42 file Click "Refinement->Refinement commands" Start "Refinement->Run refinement"

The refinement converges in 5 cycles to R(obs) 7.28% (3166 obs reflections) and wR2(all) 16.40% (8466 reflections). The final state of the refinement is available in the tutorial folder: JANA/2-dynamic/B18H22_3-refine_ori

6. Final structure

Click on "Distance calculation->Calculate distances", "Yes" to show a list of interatomic distances.

In comparison to the kinematical refinement, the standard uncertainties are now all about 0.01 Å. Bond lengths range from 1.183 to 1.451 Å.

The displacement parameters of the hydrogen atoms are now more even than with the kinematical refinement. You can see in the table below that the hydrogen atoms bonded by two electron two center bond have the ADPs about 0.055 Å but the ADPs of the hydrogen atoms bonded by two electron three center bonds have the ADPs nearly twotimes lower (0.031 Å). This difference stems from the higher draft of the bonding electrons by the boron backbone from these electrons, creating a more positive partial charge on these hydrogen atoms, which results in a higher scattering power of these hydrogens. On the other hand, the boron atoms have their ADPs about 0.0095 Å with an exception of B3, which is the boron atom surrounded by three hydrogen atoms. Its ADP value is nearly doubled in comparison to the other boron atoms, so the

effect here is opposite, which is probably due to a higher electron density around this atom, and thus its lower scattering power. You can compare these results with the kinematical refinement and with the neutron data measured at room temperature published by G.S. Pawley in Acta Crystallorg. **20** (1966) 631. Anisotropic betas were convterted to Uiso. In the neutron data, there are no similar effects on the values of the ADPs observed in the dynamically refined electron diffraction data.

When we look at the B-H distances of B3, we can see that for kinematical refinement the distances to H1b3b8 (1.21(4) Å) and H1b3 (1.16(4) Å) are within 1σ the same. The distance to H1b3b7 is 1.33(4) Å. But for the dynamical refinement, the distances are 1.450(10), 1.280(13), and 1.419(11) Å, respectively, and thus it is easy to determine, by which type of bond is the hydrogen atom to the boron attached.

Dynamical refinement						Kinematical refinement											
	Atom	Ox.	Wyck.	Site	S.O.F.	x/a	y/b	z/c	U [Å ²]	Atom	Ox.	Wyck.	Site S.O.F.	x/a	y/b	z/c	U [Å ²]
	B1		8e	1		0.43212	0.35366	0.39516	0.0082	B1		8e	1	0.43445	0.35487	0.39565	0.0033
	B2		8e	1		0.68052	0.44819	0.41197	0.0068	B2		8e	1	0.68030	0.44880	0.41293	0.0014
	B3		8e	1		0.54485	0.53529	0.20017	0.0179	B3		8e	1	0.54349	0.53238	0.20201	0.0057
	B4		8e	1		0.58894	0.33778	0.34732	0.0093	B4		8e	1	0.58898	0.33734	0.34794	0.0075
	B5		8e	1		0.65497	0.44078	0.24862	0.0092	B5		8e	1	0.65511	0.43798	0.24796	0.0084
	B6		8e	1		0.55762	0.38351	0.50110	0.0084	B6		8e	1	0.55899	0.38277	0.50211	0.0048
	B7		8e	1		0.64288	0.57084	0.32828	0.0092	B7		8e	1	0.64344	0.57126	0.32784	0.0093
	B8		8e	1		0.49517	0.39300	0.22952	0.0133	B8		8e	1	0.49569	0.39647	0.23019	0.0109
	B9		8e	1		0.57525	0.53391	0.49383	0.0065	B9		8e	1	0.57658	0.53342	0.49288	0.0059
	H1b7		8e	1		0.71869	0.65645	0.32324	0.0412	H1b3b8	3	8e	1	0.44094	0.49575	0.22316	0.0251
	H1b3b7		8e	1		0.53500	0.61710	0.29490	0.0306	H1b3		8e	1	0.54157	0.58136	0.10883	0.0281
	H1b1b9		8e	1		0.36134	0.44482	0.38948	0.0304	H1b7		8e	1	0.71616	0.64759	0.31492	0.0482
	H1b1		8e	1		0.36304	0.26575	0.41089	0.0552	H1b4		8e	1	0.62022	0.23213	0.31854	0.0220
	H1b5		8e	1		0.73941	0.41790	0.18312	0.0530	H1b2		8e	1	0.79404	0.45559	0.45447	0.0632
	H1b6		8e	1		0.58921	0.32185	0.58637	0.0635	H1b3b7	, ,	8e	1	0.53624	0.60695	0.29368	0.0221
	H1b8		8e	1		0.45452	0.33407	0.13840	0.0565	H1b8		8e	1	0.46070	0.32564	0.12515	0.0701
	H1b4		8e	1		0.61383	0.24142	0.33650	0.0561	H1b5		8e	1	0.73743	0.39968	0.18340	0.0283
	H1b2		8e	1		0.78679	0.44754	0.45536	0.0511	H1b1b9)	8e	1	0.36127	0.44474	0.39803	0.0187
	H1b3b8		8e	1		0.42448	0.48890	0.23075	0.0313	H1b6		8e	1	0.58648	0.32039	0.59949	0.0336
	H1b3		8e	1		0.53647	0.58761	0.09678	0.0599	H1b1		8e	1	0.36606	0.26634	0.41287	0.0480



| 🖬 🖶 | 🚼 🗢 ∧ 🗝 🛞 🍸 💷 🗶 | Ø 🔀 🍒 💠 🖽 🖓 | 👗 🞯 🔊 (?) | 🔤 🎇 ∪ | 🕾

Pawley 1966, 293 K, neutron single crystal											
Atom	Ox.	Wyck.	Site	S.O.F.	x/a	y/b	z/c	U [Å ²]			
B1		8e	1		-0.06590	0.14480	0.10390	0.0204			
B2		8e	1		0.17870	0.05320	0.08700	0.0191			
B3		8e	1		0.04580	-0.03310	0.29800	0.0244			
B4		8e	1		0.08820	0.16230	0.15130	0.0208			
B5		8e	1		0.15410	0.06170	0.25020	0.0225			
B6		8e	1		0.05800	0.11630	-0.00050	0.0182			
B7		8e	1		-0.07590	0.03180	-0.00590	0.0173			
B8		8e	1		0.14410	-0.06750	0.17180	0.0216			
B9		8e	1		-0.00300	0.10590	0.26990	0.0230			
H1b3b8		8e	1		-0.06290	0.01170	0.26380	0.0227			
H1b3		8e	1		0.03940	-0.07620	0.39030	0.0307			
H1b7		8e	1		0.21250	-0.13550	0.17680	0.0264			
H1b4		8e	1		0.12100	0.24920	0.16380	0.0247			
H1b2		8e	1		0.26930	0.05830	0.05160	0.0217			
H1b3b7		8e	1		0.04520	-0.10380	0.20920	0.0226			
H1b8		8e	1		-0.04080	0.15940	0.35190	0.0287			
H1b5		8e	1		0.23150	0.09110	0.31750	0.0271			
H1b1b9		8e	1		-0.13020	0.05040	0.09810	0.0190			
H1b6		8e	1		0.08720	0.16410	-0.07880	0.0208			
H1b1		8e	1		-0.13060	0.21380	0.08800	0.0252			